

# Decoding CODA files for the Moller polarimeter

Don Jones

Dec. 2019

# Decoding the binary CODA files

This was accomplished by a combination of

1. Using cefdmp on a Hall A computer to view CODA moller\_data\_XXXXX.dat files event by event
2. Reading the binary in hex mode in a viewing program (Emacs M-x hexl-mode)
3. Running the standalone CODA event viewer code provided by Bob Michaels  
[http://hallaweb.jlab.org/software/tools/simpleAna\\_18Oct2016.tar](http://hallaweb.jlab.org/software/tools/simpleAna_18Oct2016.tar)
4. Finding snippets of explanations about CODA online such as  
[https://hallaweb.jlab.org/equipment/daq/dstruct\\_year2001.html](https://hallaweb.jlab.org/equipment/daq/dstruct_year2001.html)
5. Tracing the logic/action of the FORTRAN analyzer rawread.f
6. User manuals for the TDC and scaler (included in Git repos)

# CODA header-- 32 bytes (8 words)

NOTE!! A CODA header event starts every block=0x8000 bytes!

```
0020 0000 0000 0000 0800 0000 0800 0000  
0020 0000 0200 0000 aa00 0000 0001 dac0
```

c0da0100 in little endian

## Header word definitions

1. Size of block in 4-byte words (=8192)
2. Block number starting at 0
3. Size of header in 4-byte words (=8)
4. Start of first event in this block
5. Number of words used in this block
6. Version of file format (=1)
7. Reserved
8. Magic number for error detection (=c0da0100)

Note: these headers are excluded from cefdmp output

# CODA data types and event types

word\_1 gives length of event in  
4-byte words (not including itself)

0400 0000 cc01 1100 1a31 4f5d 7443 0000  
0100 0000

word\_2 encodes event type and data type. If the  
word = aabb ccdd (which is little endian for  
0xddccbbaa)

- event type = word >> 16 (or ddcc)\*\*
- data type = (word & 0x0000ff00) >> 8 (or bb)


In the example above we have an event type of  
0x11=17 and a data type of 0x1=1

\*\*dd always seems to be 0

Name	Data type	Event type	Word 2	Word 2 in Hex
<b>Sync</b>	1	16	cc01 1000	0x1001cc
<b>Prestart</b>	1	17	cc01 1100	0x1101cc
<b>Go</b>	1	18	cc01 1200	0x1201cc
<b>Pause</b>	1	19	cc01 1300	0x1301cc
<b>End</b>	1	20	cc01 1400	0x1401cc
<b>EPICS</b>	16	131	cc10 8300	0x8310cc
<b>Physics++</b>	16	<15	cc10 0100	0x110cc


++Event types <15 are various types of triggers

# Decoding CODA data type=1, event type=16 (Sync)




Note! This is **not an actual sync event** example. I didn't come across a "sync" event in the files I looked at but inferred from the kumac that the only interesting information from a sync event is the current time directly following the type word.

# Decoding CODA data type=1, event type=17 (Prestart)



Actual prestart event example.


# Decoding CODA data type=1, event type=18 (Go)



Actual Go event example.

The current time from the Go event is the closest we have to start of data taking time.


# Decoding CODA data type=1, event type=19 (Pause)



Note! This is **not an actual pause** event example. I didn't come across a pause event in the files I looked at but inferred from the kumac that the only interesting information from a pause event is the current time directly following the type word.



# Decoding CODA data type=1, event type=20 (End)



Actual End event example.

Note that there are still events that occur after the End event hits the datastream. **However, events after the End event are not utilized.**

# Decoding CODA data type=16, event type=131 (EPICS)

Beginning of EPICS event example.


```
e802 0000 cc10 8300 e602 0000 0003 0000
5361 7420 4175 6720 3130 2031 373a 3033 or Sat Aug 10 17:03
3a35 3920 4544 5420 3230 3139 0a49 504d or :59 EDT 2019.IPM
3143 3230 2e58 504f 5320 2020 2020 2020 or 1C20.XPOS
2020 2020 2020 2020 2020 202d 302e 3033 or -0.03
```

EPICS events are recorded in ASCII not binary.

No one has built a decoder for Moller EPICS events yet to my knowledge.

# Decoding CODA data type=16, event type=1 (Physics)

## ADC readout



### cefdmp terminal output for this event

```
0x1 length=47 bank
0xc000 length=5 integer
1 1 0
0x7 length=40 integer
36 24 32 56
405 56 56 52
389 21 409 385
739 385 22 15
31 20 15 38
22 24 15 12
17 22 5 65716
196785 327865 720932 1507483
5 61665 720898 1889608
65408 7187
```

\*\*ROC number = (0x70101 & 0xff0000)>>16=7

# Specific ADC event example

Note: the meaning of words marked with '?' is not known but these are not utilized in the analysis chain.

event len	CODA type	headr len	?
event num	?	evnt stat	ROC len
ROC Num	Moll type	N ADC cha	ADC [0]
ADC [1]	ADC [2]	ADC [3]	ADC [4]
ADC [..]	ADC[..]	ADC [N-1]	N TDC cha
TDC [0]	TDC [1]	TDC [2]	
TDC[N-1]	stat len	trig info	?
itick	?	rtrn code	

2e00 0000	cc10 0100	0400 0000	0001 00c0
0100 0000	0100 0000	0000 0000	2700 0000
0101 0700	2400 0000	1800 0000	2000 0000
3800 0000	9501 0000	3800 0000	3800 0000
3400 0000	8501 0000	1500 0000	9901 0000
8101 0000	e302 0000	8101 0000	1600 0000
0f00 0000	1f00 0000	1400 0000	0f00 0000
2600 0000	1600 0000	1800 0000	0f00 0000
0c00 0000	1100 0000	1600 0000	0500 0000
b400 0100	b100 0300	b900 0500	2400 0b00
9b00 1700	0500 0000	e1f0 0000	0200 0b00
48d5 1c00	80ff 0000	131c 0000	

# General ADC event structure

Evt Len (including header)	Evt type	Header Len=-4	?
Evt. num.	?	Evt status	ROC Len
Coda info (including ROC)	Moller evt type	N_ADC ch	ADC[0]
...	...	ADC[N_ADC-1]	N_TDC ch
TDC[0]	...	TDC[N_TDC-1]	N_TDC2 ch**
TDC2[32]	...	TDC2[N_TDC2]	Status Record Len=5
Trigger info (itrig)	?	itick	?
iret			

\*\*This second bank of TDC information is only read out if Physics event type=37. I don't have an example of this type of event.

TDC information: (see Scaler.pdf in Git repos)

- The first word is N\_TDC
- The TDC is only 16 bit which only requires half the 32 bit word so the first 16 bits are the TDC values TDC[i] & 0xffff which for the example below gives 0xb4, 0xb1, 0xb9, 0x24, 0x9b
- The last 17th bit encodes the TDC trigger phase (ited) with leading edge = 1 and trailing edge = 0
  - $ited[i] = (TDC[i] \gg 16) \& 1$ , which for the first example below is  $(0x100b4 \gg 16) \& 1 = 1$
- Bits 18-22 are channel number (itch)
  - $itch[i] = ((TDC[i] \gg 17) \& 0x1f)$  which for the example below is 0, 1, 2, 5, 11 (Note that the FORTRAN code increments these by 1 so the channel numbering starts at 1)

```
0500 0000 b400 0100 b100 0300 b900 0500
2400 0b00 9b00 1700
```


Trigger information:

- The first word of the status record is the length of the status record
- The first 8 bits of the second word is the 8 channels of itrig (either 0 or 1) encoded bitwise so 0xe1 = 1110 0001 giving itrig[]={1,0,0,0,0,1,1,1}
- itick is the optional 4th word of the status record and is the number of ticks per event of the 120 Hz CPU clock
- iret is the optional 6th word of the status record and is a return/error code

```
0500 0000 e1f0 0000 0200 0b00 48d5 1c00  
80ff 0000 131c 0000
```

# Decoding CODA data type=16, event type=1 (Physics)

## Scaler readout





trig info  
?

## Specific scaler readout example

3300	0000	cc10	0100	0400	0000	0001	00c0
0405	0000	0100	0000	0000	0000	2c00	0000
0401	0700	2000	0000	2000	0000	130e	0100
fdbd	0000	2b2e	0000	1900	0000	523c	0000
222e	0000	0000	0000	0000	0000	0000	0000
d3ff	0100	f63e	0000	d7ff	0300	0000	0000
0000	0000	d9ff	0300	0000	0000	130e	0100
fdbd	0000	2b2e	0000	1800	0000	110e	0100
fdbd	0000	f6d3	0200	0000	0000	8202	0000
1709	0a00	a6a0	0a00	1b96	0a00	a95d	0a00
0000	0000	0000	0000	8202	0000	2400	0000
0000	0000	0000	0000	0500	0000	e2f0	0000
0000	0100	89d6	1c00	84ff	0000	0000	0000

event len	CODA type	headr len	?
event num	?	evnt stat	ROC len
ROC Num	Moll type	N Scal ch	Sca [0]
Sca [0]	Sca [0]	Sca [0]	Sca [0]
Sca [..]	Sca [..]	Sca [N-1]	Moll type
N_ADC cha	N_TDC_cha	stat len	trig info
?	itick	?	rtrn code

```
**ROC = (0x70104 & 0xff0000)>>16=7
```

# General scaler event structure

Evt Len (including header)	Evt type	Header Len=4	?
Evt. num.	?	Evt status	ROC Len
ROC info	Moller evt type	N_Sca ch=32	Sca[0]
...	...	Sca[N_Sca-1]	Moller evt type=36
N_ADC=0	N_TDC=0	Status Record Len=5	Trigger info (itrig)
?	itick	?	iret

It looks like in order to include the status record for scaler events, they append an ADC+TDC+status record event with 0 channels for ADC and TDC